



BrainScan™ Installation Details Form (IDF)

We are unable to enter purchase orders into the system or guarantee delivery dates without a technically accurate and fully completed IDF.

BrainScan 1: Includes remote Hot Water inlet, Cold inlet/Recirculated Return Water and Mixed Water outlet temperature readings (integral to DRV80). Also gives the ability to remotely change Mixed water outlet temperature setpoint. Includes systm graphic file.

BrainScan 2*: Supplied as BrainScan 1 with Hot Water inlet, Cold Water inlet and Mixed Water outlet pressure transmitters.

BrainScan 3*: Provided as BrainScan 2 with Mixed Water outlet and Re-Circulation Return flow meters.

Brain Model Number:(DRV80, 80R, DMC etc) _____

BrainScan Model: BrainScan 1 _____ BrainScan 2* _____ BrainScan 3* _____ (check one)

*Only available with Digital Mixing Centers (DMC)

System Interface Options: (Please Select One)

Building Automation System

Protocol Options: 1 - 4

- 1) LonWorks 2) Modbus RTU 3) BACnet Metasys N2 4) BACnet MSTP

For selections 2 - 4 above please provide:

Baud Rate _____ MAC address _____ Node ID _____

Protocol Options: 5 - 6

- 5) Modbus TCP 6) BACnet IP

For selections 5 - 6 above please provide:

IF Static IP address _____ Subnet Mask _____ Gateway _____ Node ID _____
 OR
 IF DHCP Node ID _____

Web Browser Interface

TCP/IP Configuration

TCP/IP address _____

Subnet Mask _____

Gateway _____

Primary DNS _____

Secondary DNS _____

DNS Domain Suffix(eg: mycompany.com) _____

E-Mail Notification Factory Enable Y N

Mail Server Name/Address _____

Mail Server Port (default 25) _____

"From" E-Mail Address _____

Notification List for Alerts

1 _____

2 _____

3 _____